

# ECG Module

Zug 3/5 leads ECG module, referred herein after as MECG35 Module is designed to be integrated in adult, children or neonate patient monitoring systems in all hospital departments and emergency services.

The MECG35 ECG Module can measure ECG, Respiratory, and Temperature.

It allows the detection of 18 types of arrhythmias.

## ECG Module

MECG35



### FEATURES

- Support 3/5 leads ECG (3 independent channels), 1 respiration, 2 temperature
- With 4 measurement mode: diagnose, monitor, HARDEST and operation measurement
- Arrhythmia diagnostic (18 types)
- With special function of ECG analysis, verified by MIT database
- Support ZUGMED protocol
- Support OEM, the best choice of multi-parameter monitor

### APPLICATION

- ECG function of Patient Monitor

## OVERVIEW

Zug 3/5 leads ECG module, referred herein after as MEG35 Module is designed to be integrated in adult, children or neonate patient monitoring systems in all hospital departments and emergency services.

The MEG35 ECG Module can measure ECG, Respiratory, and Temperature.

### Methods of measurements:

**The ECG measurement:** as the heart generates electrical stimulation and biological electricity, these signals

transmit to body surface before shrinking mechanically. Parts of the body generate different electrical potential and form the potential difference on body surface. The recorded potential difference creates dynamic curve which is ECG signal. **The respiratory signal measurement** is based on the impedance method and detects vital parameters by measuring the change of thoracic impedance when people breathe and the thorax fluctuation. **Temperature measurement** is performed through thermostat.

## FUNCTIONS

### MEASUREMENT MODES

There are four kinds of measurement mode: diagnose, monitor, HARDEST and operation mode.

**Diagnose mode** ..... The filter range is 0.05HZ-130HZ

**Monitor mode** ..... The filter range is 0.05HZ-40HZ

**HARDEST mode** ..... The filter range is 5HZ-20HZ

**Operation mode** ..... The filter range is 1HZ-25HZ

### PATIENT TYPES

Adult, Children, Neonate.

### CALIBRATION

Recording the amplitude of wave by inputting a standard voltage 1mv.

### GAIN SETTING

Adjusting the amplitude of the ECG and respiratory waveform.

### HEART RATE CALCULATION CHANNEL

I, II, VI channel can be used as analysis channel or you can choose any other channel by yourself.

### LEADS SETTING

You can select Lead I, Lead II, Lead AVL, Lead AVR, Lead AVF as the signal of Channel I and Channel II.

### RELIABILITY

The trapped wave mode: 50HZ, 60HZ, 50/60HZ are available and you also can turn off the trapped wave mode.

### ARRHYTHMIA CALCULATION LIST

| No | Description             | No | Description                  |
|----|-------------------------|----|------------------------------|
| 1  | Asystole                | 10 | Supraventricular tachycardia |
| 2  | VF                      | 11 | Supraventricular bradycardia |
| 3  | Ventricular tachycardia | 12 | Multiform VPB                |
| 4  | R ON T                  | 13 | Pace not capture             |
| 5  | Multiple VPB            | 14 | Pacer not paced              |
| 6  | Couple VPB              | 15 | Irregular rhythm             |
| 7  | Accidental VPB          | 16 | Missed beat                  |
| 8  | Bigeminy                | 17 | Normal sinus rhythm          |
| 9  | Trigeminy               | 18 | Oversized noise              |

## PERFORMANCE

### ECG

**Range** ..... 0.15mv-5.5mv

**Accuracy** ..... 2.36uV/LSB

**Resolution** ..... 1%

### RESPIRATORY

**Range** ..... 15-120 rpm

**Base resistance value** ..... 500-2000  $\Omega$

**Variable resistance value** ..... 0.2-3.0  $\Omega$

### HEART RATE

**Range** ..... 15-300 bpm, 15-350 bpm for neonate

**Accuracy** .....  $\pm 1$  bpm

**Resolution** ..... 1 bpm

### TEMPERATURE

**Range** ..... 0-50°C

**Accuracy** ..... /

**Resolution** ..... 0.1°C

## STANDARDS COMPLIANCE DESCRIPTION

| Standard NO.  | Number of standard  | Version |
|---------------|---|---------|
| IEC60601-2-25 | Medical electrical equipment—part 2-25:<br>particular requirements for the basic safety and<br>essential performance of electrocardiographs | 2011    |

## ELECTRICAL CHARACTERISTICS


|                                |  |
|--------------------------------|--|
| <b>Input Voltage</b> .....     | External power supply should provide +12V DC<br>Voltage offset range should between $\pm 10\%$ of voltage full range |
| <b>Power Consumption</b> ..... | $\leq 5W$  |

## ENVIRONMENT CHARACTERISTICS

|                    | Operating environment        | Storage environment           |
|--------------------|------------------------------|-------------------------------|
| <b>Temperature</b> | 10°C to 40°C (50°F to 140°F) | -20°C to +70°C (4°F to 158°F) |
| <b>Humidity</b>    | 15%-90% non-condensing       | 15%-90% non-condensing        |
| <b>Altitude</b>    | -170 m to +1700 m            | -170 m to +1700 m             |

## DIMENSIONS


The board dimensions are given in millimeters.



# INTERFACES

## POWER AND COMMUNICATION INTERFACE

The connector J4 shown below is used for the purpose of communication and power supply. The pin 1 is actually indicated by an arrow on the PCB.



The table shows the pin assignment of this connector.

| Pin No | Signal | Description                                 |
|--------|--------|---|
| 1      | TXD    | UART Sending data from module to host       |
| 2      | RXD    | UART Receiving data from host to the module |
| 3      | DGND   | Ground                                      |
| 4      | +12V   | Power Supply input 12V DC                   |
| 5      | DGND   | Ground                                      |
| 6      | +12V   | Power Supply input 12V DC                   |


Note: Pins 4 and 6 are internally connected to +12V, and pins 3 and 5 are internally connected to DGND.

The table shows the pin assignment of this connector.

| Pin No | Signal     | Description                               |
|--------|------------|---|
| 1      | RA         | Right arm                                 |
| 2      | LA         | Left arm                                  |
| 3      | LL         | Left leg                                  |
| 4      | V1         | Chest lead, in No.4 ribs stick to sternum |
| 5      | ECG SHIELD | Signal shielding wire                     |
| 6      | RL         | Right leg                                 |

## TEMPERATURE SENSOR INTERFACE

The connector J8 shown below is used for the purpose of temperature signals transmission. The pin 1 is actually indicated by an arrow on the PCB.




The table shows the pin assignment of this connector.

| Pin No | Signal | Description                          |
|--------|--------|--------------------------------------|
| 1      | TEMP1  | Temperature sensor 1, positive input |
| 2      | TGND   | Temperature sensor 1, negative input |
| 3      | TEMP2  | Temperature sensor 2, positive input |
| 4      | TGND   | Temperature sensor 2, negative input |

## ECG LEADS INPUT INTERFACE

The connector J10 shown below is used for the purpose of ECG signals transmission. The pin 1 is actually indicated by an arrow on the PCB.



# ORDERING

Our 3/5 Leads module part number is **MECG35**.

For ordering our module, please contact directly our sales team by email at [sales@zugmed.com](mailto:sales@zugmed.com) or refer to our website [www.zugmed.com](http://www.zugmed.com) for further information.

Headquarter  
**Zug Medical Systems SAS**  
 2000 Route Des Lucioles, CS 90029  
 06901 Sophia Antipolis  
 France

W [www.zugmed.com](http://www.zugmed.com)  
 E [sales@zugmed.com](mailto:sales@zugmed.com)  
 T +33 (0)9 78 234 726  
 F +33 (0)9 57 135 843

**ZUG**  
 MEDICAL SYSTEMS